



Advanced Bionics



A Sonova brand

TECHNICAL SPECIFICATIONS


HiRes™ Ultra 3D Cochlear Implant
with the HiFocus SlimJ and Mid-Scala Electrodes

The Foundation of Better Hearing

The new HiRes Ultra 3D Cochlear Implant is designed to provide all the benefits of the HiRes Ultra implant platform but with a hassle free MRI solution. With a patented magnet design¹ that allows the magnet to stay in place for high-resolution 3T MR scans, your patient can have an MRI procedure without any pain, and enjoy uninterrupted hearing.

Developed with leading cochlear implant surgeons and using state-of-the-art manufacturing processes, the HiRes Ultra 3D with the HiFocus™ SlimJ and HiFocus™ Mid-Scala electrodes offer unique features designed to suit individual patient anatomy and surgical preferences for the best possible hearing outcomes.

HiRes Ultra 3D Cochlear Implant



IMPLANT TECHNICAL SPECIFICATIONS		
Information Update Rate	90 kHz	
Stimulation Rate	Up to 83,000 pps (software limited)	
Independent Output Circuits	16	
Spectral Bands	Up to 120 sites of stimulation (software limited)	
Communication Link	Bi-directional inductive link	
IntelliLink™ Safety Feature	Implant and processor association	
Diagnostics	Neural response imaging (NRI), impedance measurements, ESRT, Integrity testing	
ADC Resolution Sampling Rate	Resolution: 9 bits, Sampling Rate: 25 kHz	
Pulse Amplitude & Width	0 - 2040 µA & 10.78 - 229 µS	
Impedance Accuracy	2.5 kΩ	
STIMULATION DELIVERY SPECIFICATIONS		
HiRes Fidelity 120™ Strategies	HiRes Optima*-P (paired), HiRes Optima*-S (sequential), ClearVoice™*, HiRes-P with Fidelity 120™* (paired), HiRes-S with Fidelity 120™* (sequential)	
HiResolution™ Sound Strategies	HiRes-P (paired) and HiRes-S (sequential)	
Conventional Strategies	CIS and MPS Modes	
IMPLANT MATERIALS AND DIMENSIONS		
Titanium Case	4.5 mm titanium case including silicone	
Antenna Coil	3.6 mm antenna coil silicone	
Housing	25/28.5 mm x 56.2 mm flexible silicone	
Weight	11 grams	
Volume	4800 mm ³	
Magnet	Neodymium inside a titanium case	
Telemetry Coil	Gold-braided wire and platinum-shield wire in flexible silicone, reinforced with high-density polymer fiber	
Ground	2 — Case ground and ring electrode ground	
Impact Resistance Value	Exceeds the impact requirements specified in EN45502-2-3:2010	
Pressurized Environment Information	Can withstand a pressure up to a depth of 42m under water (138 feet) or a gauge pressure of 4ATM (413kPa)	
ELECTRODE TECHNICAL SPECIFICATIONS		
	HiFocus SlimJ	HiFocus Mid-Scala
Electrodes	16 platinum contacts; platinum-iridium wires; flexible silicone carrier; integrated ground on lead	
Minimum Exposed Contact Area	0.12 mm ²	0.12 mm ²
Contact Spacing	1.3 mm	0.975 mm
Active Length	~20 mm	~15 mm
Forceps	Yes	Jeweler's forceps or a similar instrument can be used
Insertion Tool	–	Optional
Freehand	Yes	Yes
Reloadable	–	2 Maximum
Insertions	3 Maximum	3 Maximum
Recommended Insertion	Round window, extended round window and cochleostomy	
SURGICAL PROCEDURES – SEE SYSTEM IFU (REF) FOR DETAILS		
Magnet	Can be left in place for 3T and 1.5T MRI scans (no bandaging required) or can be easily removed for 3T or 1.5T MRI scans. See the MRI Safety Information booklet.	
Electrosurgical Instruments	Monopolar cautery can be used outside of the head and neck, bipolar 1mm from the implant.	
Extracorporeal Electrical Stimulation	Probes not over the implant.	
Neurostimulation	Not directly over the implant.	
Diathermy	Safe to use shortwave and longwave.	
Radiation	250 Gray using a 15MeV beam strength and 3cm depth. Diagnostic radiation imaging such as CT, X-Ray, Mammography, etc. are safe to use.	
Diagnostic Ultrasound	Safe to use.	

References

1. Lee et al. (2018). *U.S. Patent No. 9,919,152B2*. Valencia, CA: U.S. Patent and Trademark Office.

*Not approved for pediatric use in the United States.



Advanced Bionics AG

Laubisrütistrasse 28,
8712 Stäfa, Switzerland
T: +41.58.928.78.00
F: +41.58.928.78.90
info.switzerland@AdvancedBionics.com

Advanced Bionics LLC

28515 Westinghouse Place
Valencia, CA 91355, United States
T: +1.877.829.0026
T: +1.661.362.1400
F: +1.661.362.1500
info.us@AdvancedBionics.com

For information on additional AB locations, please visit
[AdvancedBionics.com/contact](https://www.AdvancedBionics.com/contact)

AB – A Sonova brand

Please contact your local AB representative for
regulatory approval and availability in your region.